

# Priming Volumes - approximate

<i>Double Lumen</i>	PROXIMAL	DISTAL	<i>Triple Lumen</i>	PROXIMAL	MEDIAL	DISTAL
<b>4FR x 8cm</b>	0.25	0.23	<b>5FR x 8cm</b>	0.25	0.25	0.25
<b>4FR x 13cm</b>	0.25	0.23	<b>5FR x 13cm</b>	0.30	0.25	0.25
<b>4FR x 30cm</b>	0.32	0.30	<b>5FR x 30cm</b>	0.35	0.32	0.40
<b>5FR x 8cm</b>	0.25	0.22	<b>7FR x 15cm</b>	0.40	0.36	0.43
<b>5FR x 13cm</b>	0.27	0.25	<b>7FR x 20cm</b>	0.40	0.38	0.47
<b>5FR x 25cm</b>	0.35	0.32	<b>7FR x 30cm</b>	0.45	0.40	0.57
<b>7FR x 15cm</b>	0.41	0.45				
<b>7FR x 20cm</b>	0.45	0.52				
<b>7FR x 30cm</b>	0.50	0.66	<b>QUAD -SA8FQ20</b>	18ga White 0.38ml	18ga Blue 0.37ml 14ga Gray 0.46ml	16ga Brown 0.41ml